

**GRO $\gamma$ /MGSA (C) Anti-Human Rabbit IgG Affinity Purify**

**Polyclonal Antibody**

**Cat.NO.: PA01188**

---

3th Edition

**Description:** GRO  $\gamma$ /MGSA as IL-8 family has the high homology with IL-8 in amino acid sequences. It is reported that it also has the similar configuration in 3D structural analysis and the almost same level of leukocyte chemotaxis with IL-8.

**Antigen:** Synthetic peptide in portion of C-terminus of Human GRO $\gamma$ /MGSA.

**Form:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1 $\mu$ g/ml.

**Specificity:** Non-Crossreacts to Human GRO $\gamma$ , Human IL-8, Human IP-10, Rat GRO/CINC-1, Rat GRO/CINC-2 $\gamma$ , Rat GRO/CINC-2 $\beta$  and Rat GRO/CINC-3.