

GRO α (C) Anti-Human Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01189

3th Edition

Description:GRO (Growth Related Oncogene) belonging to IL-8 family is polypeptide which has three isoforms α , β and γ , and inhibits proliferation of endothelial cells. It is reported that only GRO α inhibits angiogenesis in case of administration of GRO in vivo.

Antigen:Synthetic peptide in portion of C-terminus of Human GRO α

Form:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1 μ g/ml.

Specificity:Recombinant Human GRO α Specific. Non-crossreacts Recombinant Human GRO β