

Syndecan-4 Anti-Human Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01196

3th Edition

Description: Syndecan, a cell surface proteoglycan, is an integral membrane protein acting as a receptor for the extracellular matrix. The syndecan family of four transmembrane heparan sulfate proteoglycans binds a variety of soluble and insoluble extracellular effectors.

Antigen: Maltose Binding Protein (MBP) - Syndecan-4 recombinant protein

Form: Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Confirmed by western blotting