

Thioredoxin (C) Anti-Human Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01201

3th Edition

Description: This antibody can be used for western blotting in concentration of 2 ~ 5 µg/ml.

Antigen: Synthetic peptide in portion of C-terminus of Human Thioredoxin

Form: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1-5 µg/ml.

Specificity: