

Anti-Human IgG (Fc) Goat IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01256

3th Edition

Description:Anti-Human IgG (Fc) Goat IgG Affinity Purify, 50 µg. Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃, purified with antigen peptide.

Antigen:Human IgG

Form:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃.

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:The antibody is absorbed by goat IgG, rabbit IgG, mouse IgG and rat IgG.