

**IL-1 $\beta$  (C) Anti-Rat Rabbit IgG Affinity Purify**

**Polyclonal Antibody**

**Cat.NO.: PA01260**

---

3th Edition

**Description:** The Interleukin-1 (IL-1) species represent an important family of biologically active mono nuclear cell-derived proteins which are involved in inflammatory reactions and in immune responses.

**Antigen:** Synthetic peptide to C-terminal of Rat IL-1 $\beta$

**Form:** Lyophilized product from 1% BSA in PBS containing 0.05%NaN<sub>3</sub>

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1 $\mu$ g/ml.

**Specificity:**