

IL-1 β (C) Anti-Rat Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01260

3th Edition

Description: The Interleukin-1 (IL-1) species represent an important family of biologically active mono nuclear cell-derived proteins which are involved in inflammatory reactions and in immune responses.

Antigen: Synthetic peptide to C-terminal of Rat IL-1 β

Form: Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1 μ g/ml.

Specificity: