

**IL-6 (N) Anti-Rat Rabbit IgG Affinity Purify**

**Polyclonal Antibody**

**Cat.NO.: PA01262**

---

3th Edition

**Description:** Interleukin-6 (IL-6) is one of cytokines which was found as a B cell differentiation-inducing factor. In these days, IL-6 is generating a lot of attention not only from antibody-producing system but also from assorted physiologically activity such as induction in biosynthesis of acute stage protein in hepatocytes and acceleration of proliferation based on synergism with IL-3.

**Antigen:** Synthetic peptide in portion of N-terminus of Rat IL-6

**Form:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**