

Anti-Rab27A/B Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01305

3th Edition

Description: Small GTPase Rab is a large family of membrane trafficking proteins that are conserved in all eukaryotic cells. More than 60 Rab isoforms have been reported in mice and humans, and they are believed to regulate various steps (or various types) of organelle transport.

Antigen: Recombinant Rab27B protein

Form: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Cross-react with mouse and rat Rab27A/B. Not cross-react with the Rab3 subfamily