

Rictor (R1291) Anti-Human Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01330

3th Edition

Description:Rictor (rapamycin-insensitive companion of mTOR/ mAVO3) is about 200 kDa and identified as an mTOR-binding protein. Rictor is a factor that constitutes a complex with mLST8 and Sin 1 which is called mTORC2 (mTOR complex 2), and it is contrasted with raptor in mTORC1.

Antigen:Synthetic peptide of a part of Human Rictor (GSSHTLPRRAQSLKA)

Form:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: