

**Nephrin (N) Anti-Human Rabbit IgG Affinity Purify**

**Polyclonal Antibody**

**Cat.NO.: PA01337**

---

3th Edition

**Description:** Nephrin has been identified as a disease causing molecule responsible for congenital nephrotic syndrome of Finnish type, in which the proteinuria starts in utero or immediately after birth (ref. 1).

**Antigen:** Recombinant protein of Human Nephrin

**Form:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Please see datasheet.