

Septin 7 (C) Anti-Human Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01361

3th Edition

Description: A cytoskeletal or scaffold protein family named “septin” has been shown to be involved in diverse phenomena in the living body, including cell division. The term “septin” is derived from “septation,” which means division.

Antigen: Synthetic peptide of the C terminal part of Human Septin 7 (NSSRTLEKNKKKKGKIF)

Form: Lyophilized product from 1 % BSA in PBS containing 0.05% NaN₃

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Cross reacts with mouse and rat.