

Anti-Human/Mouse/Rat Phospho-CANX (Ser583) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01516

3th Edition

Description: Calcium-binding protein that interacts with newly synthesized glycoproteins in the endoplasmic reticulum. It may act in assisting protein assembly and/or in the retention within the ER of unassembled protein subunits. It seems to play a major role in the quality control apparatus of the ER by the retention of incorrectly folded proteins.

Antigen: Synthesized peptide derived from human Calnexin around the phosphorylation site of Ser583

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: