

Anti-Human/Mouse/Rat Phospho-CASP6 (Ser257) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01517

3th Edition

Description:Involved in the activation cascade of caspases responsible for apoptosis execution. Cleaves poly(ADP-ribose) polymerase in vitro, as well as lamins. Overexpression promotes programmed cell death.

Antigen:Synthesized peptide derived from human Caspase-6 around the phosphorylation site of S257.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: