

Anti-Human Phospho-CDC25C (Ser216) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01523

3th Edition

Description: Functions as a dosage-dependent inducer in mitotic control. It is a tyrosine protein phosphatase required for progression of the cell cycle. It directly dephosphorylates CDK1 and activate its kinase activity.

Antigen: Synthesized peptide derived from human Cdc25C around the phosphorylation site of Ser216

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: