

**Anti-Human/Mouse/Rat Phospho-c-Myc (Thr58) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01531**

---

3th Edition

**Description:** Participates in the regulation of gene transcription. Binds DNA in a non-specific manner, yet also specifically recognizes the core sequence 5'-CAC[GA]TG-3'. Seems to activate the transcription of growth-related genes.

**Antigen:** Synthesized peptide derived from human c-Myc around the phosphorylation site of Thr58

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**