

Anti-Human/Mouse/Rat Phospho-FLK1 (Tyr1175) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01555

3th Edition

Description: Receptor for VEGF or VEGFC. Has a tyrosine-protein kinase activity. The VEGF-kinase ligand/receptor signaling system plays a key role in vascular development and regulation of vascular permeability. In case of HIV-1 infection, the interaction with extracellular viral Tat protein seems to enhance angiogenesis in Kaposi's sarcoma lesions.

Antigen: Synthesized peptide derived from human Flk-1 around the phosphorylation site of Tyr1175

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Detected in cornea (at protein level). Widely expressed.