

Anti-Human/Mouse/Monkey Phospho-GAB1 (Tyr627) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01562

3th Edition

Description: The protein encoded by this gene is a member of the IRS1-like multisubstrate docking protein family. It is an important mediator of branching tubulogenesis and plays a central role in cellular growth response, transformation and apoptosis. Two transcript variants encoding different isoforms have been found for this gene.

Antigen: Synthesized peptide derived from human Gab 1 around the phosphorylation site of Tyr627

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: