

Anti-Human/Mouse/Rat/Monkey Phospho-IRS1 (Ser1101) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01577

3th Edition

Description: May mediate the control of various cellular processes by insulin. When phosphorylated by the insulin receptor binds specifically to various cellular proteins containing SH2 domains such as phosphatidylinositol 3-kinase p85 subunit or GRB2. Activates phosphatidylinositol 3-kinase when bound to the regulatory p85 subunit.

Antigen: Synthesized peptide derived from human IRS-1 around the phosphorylation site of Ser1101

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: