

Anti-Human/Mouse Phospho-SHB (Tyr246) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01629

3th Edition

Description: SHB (SH2 Domain Containing Adaptor Protein B) is a Protein Coding gene. Diseases associated with SHB include Sulfhemoglobinemia. Among its related pathways are Development VEGF signaling via VEGFR2 - generic cascades and RET signaling. GO annotations related to this gene include SH3/SH2 adaptor activity. An important paralog of this gene is ENSG00000255872.

Antigen: Synthesized peptide derived from human Shb around the phosphorylation site of Tyr246

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Widely expressed.