



Anti-Human/Mouse/Rat/Monkey Phospho-STAT4 (Tyr693) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01633

3th Edition

Description: Carries out a dual function: signal transduction and activation of transcription. Involved in IL12 signaling.

Antigen: Synthesized peptide derived from human Stat4 around the phosphorylation site of Tyr693

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: