

Anti-Human/Mouse/Rat Phospho-CSK (Ser364) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01656

3th Edition

Description: Specifically phosphorylates 'Tyr-504' on LCK, which acts as a negative regulatory site. Can also act on the LYN and FYN kinases.

Antigen: Synthesized peptide derived from human Csk around the phosphorylation site of Ser364

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in lung and macrophages.