

**Anti-Human/Mouse/Rat Phospho-CSK (Ser364) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01656**

---

3th Edition

**Description:** Specifically phosphorylates 'Tyr-504' on LCK, which acts as a negative regulatory site. Can also act on the LYN and FYN kinases.

**Antigen:** Synthesized peptide derived from human Csk around the phosphorylation site of Ser364

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expressed in lung and macrophages.