

Anti-Human DEFB132 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03031

3th Edition

Description: Defensins are cysteine-rich cationic polypeptides that are important in the immunologic response to invading microorganisms. The protein encoded by this gene is secreted and is a member of the beta defensin protein family. This protein binds spermatozoa and has antimicrobial activity against *E. coli*. Beta defensin genes are found in several clusters throughout the genome, with this gene mapping to a cluster at 20p13.

Antigen: Recombinant protein of human DEFB132

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: