

Anti-Human DHRS4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03526

3th Edition

Description:DHRS4 (Dehydrogenase/Reductase 4) is a Protein Coding gene. Among its related pathways are Signaling by Retinoic Acid and the visual cycle I (vertebrates). GO annotations related to this gene include receptor binding and oxidoreductase activity, acting on NAD(P)H, quinone or similar compound as acceptor. An important paralog of this gene is DHRS4L2.

Antigen:Synthesized peptide derived from the C-terminal region of human DHRS4

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Isoform 1 is predominantly expressed in normal cervix (at protein level). Isoform 4 is expressed in some neoplastic cervical tissues, but not in normal cervix (at protein level). Isoform 5 and isoform 6 are expressed in a few neoplastic cervical tissues.