

Anti-Human/Mouse/Rat FRK Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA12081

3th Edition

Description:The protein encoded by this gene belongs to the TYR family of protein kinases. This tyrosine kinase is a nuclear protein and may function during G1 and S phase of the cell cycle and suppress growth.

Antigen:Recombinant protein of human FRK

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Restricted to cells lines derived from tissues of lymphoid, brain, breast, colon and bladder origin.