

Anti-Human/Mouse/Rat FLCN Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA12085

3th Edition

Description: This gene is located within the Smith-Magenis syndrome region on chromosome 17. Mutations in this gene are associated with Birt-Hogg-Dube syndrome, which is characterized by fibrofolliculomas, renal tumors, lung cysts, and pneumothorax. Alternative splicing of this gene results in two transcript variants encoding different isoforms.

Antigen: Recombinant protein of human FLCN

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in most tissues tested, including skin, lung, kidney, heart, testis and stomach.