

**Anti-Human/Mouse EPHA10 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA12141**

---

3th Edition

**Description:** Ephrin receptors, the largest subfamily of receptor tyrosine kinases (RTKs), and their ephrin ligands are important mediators of cell-cell communication regulating cell attachment, shape, and mobility in neuronal and epithelial cells

**Antigen:** Recombinant protein of human EPHA10

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Mainly expressed in testis.