

**Anti-Human/Mouse/Rat DDIT4L Polyclonal Antibody**  
**Polyclonal Antibody**

**Cat.NO.: PA12238**

---

3th Edition

**Description:**REDD-2 (regulated in development and DNA damage response 2), also designated Rtp801L or DDIT4L (DNA-damage-inducible transcript 4-like), is a 193 amino acid cytoplasmic protein belonging to the DDIT4 family and is predominantly expressed in skeletal muscle. Considered a stress-induced protein, REDD-2 is a negative regulator of the mTOR (mammalian target of rapamycin) pathway. mTOR is a serine/threonine kinase that plays an essential role in cell growth control and is an important regulator of skeletal muscle size. Highly expressed in human atherosclerotic lesions and macrophages, REDD-2 mediates monocyte cell death through reduction of Trx (thioredoxin-1) expression. REDD2 expression in macrophages increases oxidized LDL (oxLDL)-induced cell death, suggesting that REDD2 may play a critical role in arterial pathology.

**Antigen:**Recombinant protein of human DDIT4L

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Up-regulated in atherosclerotic plaques relative to healthy segments of the same artery.