

**Anti-Human/Mouse/Rat HSP75 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA12319**

---

3th Edition

**Description:**HSP90 proteins are highly conserved molecular chaperones that have key roles in signal transduction, protein folding, protein degradation, and morphologic evolution. HSP90 proteins normally associate with other cochaperones and play important roles in folding newly synthesized proteins or stabilizing and refolding denatured proteins after stress. TRAP1 is a mitochondrial HSP90 protein. Other HSP90 proteins are found in cytosol and endoplasmic reticulum

**Antigen:**Recombinant protein of human TRAP1

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Found in skeletal muscle, liver, heart, brain, kidney, pancreas, lung and placenta.