

**Anti-Human/Mouse ING2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA12352**

---

3th Edition

**Description:** This gene is a member of the inhibitor of growth (ING) family. Members of the ING family associate with and modulate the activity of histone acetyltransferase (HAT) and histone deacetylase (HDAC) complexes and function in DNA repair and apoptosis

**Antigen:** Recombinant protein of human ING2

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Widely expressed. Higher expressed in colon-cancer tumor than in normal colon tissues.