

Anti-Human/Mouse Oct-1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA12526

3th Edition

Description:The OCT1 transcription factor was among the first identified members of the POU transcription factor family. Members of this family contain the POU domain, a 160-amino acid region necessary for DNA binding to the octameric sequence ATGCAAAT.

Antigen:Recombinant protein of human POU2F1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Ubiquitous. Isoform 2 is lymphocyte-specific.