

Anti-Human/Mouse DCT Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA12827

3th Edition

Description: This protein can be involved in regulating eumelanin and pheomelanin levels. The subcellular location of this protein is melanosome membrane, Single-pass type I membrane protein. It may form a multienzyme complex together with Tyrosinase and TYRP1.

Antigen: Recombinant protein of human DCT

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: