

**Anti-Human/Mouse DCT Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA12827**

---

3th Edition

**Description:** This protein can be involved in regulating eumelanin and pheomelanin levels. The subcellular location of this protein is melanosome membrane, Single-pass type I membrane protein. It may form a multienzyme complex together with Tyrosinase and TYRP1.

**Antigen:** Recombinant protein of human DCT

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**