

**Anti-Human SCGB1A1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA12846**

---

3th Edition

**Description:** This gene encodes a member of the secretoglobin family of small secreted proteins. The encoded protein has been implicated in numerous functions including anti-inflammation, inhibition of phospholipase A2 and the sequestering of hydrophobic ligands. Defects in this gene are associated with a susceptibility to asthma.

**Antigen:** Recombinant protein of human SCGB1A1

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Clara cells (nonciliated cells of the surface epithelium of the pulmonary airways).