

**Anti-Human/Mouse/Rat STXBP2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA12850**

---

3th Edition

**Description:** This gene encodes a member of the STXBP/unc-18/SEC1 family. The encoded protein is involved in intracellular trafficking, control of SNARE (soluble NSF attachment protein receptor) complex assembly, and the release of cytotoxic granules by natural killer cells. Mutations in this gene are associated with familial hemophagocytic lymphohistiocytosis. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

**Antigen:** Recombinant protein of human STXBP2

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Placenta, lung, liver, kidney and pancreas, as well as in peripheral blood lymphocytes.