

Anti-Human/Mouse/Rat AASDHPPT Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA12865

3th Edition

Description: The protein encoded by this gene is similar to *Saccharomyces cerevisiae* LYS5, which is required for the activation of the alpha-amino adipate dehydrogenase in the biosynthetic pathway of lysine. Yeast alpha-amino adipate dehydrogenase converts alpha-biosynthetic-amino adipate semialdehyde to alpha-amino adipate. It has been suggested that defects in the human gene result in pipecolic acidemia.

Antigen: Recombinant protein of human AASDHPPT

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Detected in heart, skeletal muscle, placenta, testis, brain, pancreas, liver and kidney.