

Anti-Human ACTL8 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA12897

3th Edition

Description: This protein belongs to the actin family. Actin-like 8 is a protein in humans that is encoded by the ACTL8 gene. The high levels expression in testis and pancreas and weak expression in placenta. The subcellular location of this protein in cytoplasm, cytoskeleton.

Antigen: Recombinant protein of human ACTL8

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Strongly expressed in testis and pancreas. Weak expression in placenta.