

**Anti-Human/Mouse ACTN3 Polyclonal Antibody**  
**Polyclonal Antibody**

**Cat.NO.: PA12898**

---

3th Edition

**Description:** This gene encodes a member of the alpha-actin binding protein gene family. The encoded protein is primarily expressed in skeletal muscle and functions as a structural component of sarcomeric Z line. This protein is involved in crosslinking actin containing thin filaments. An allelic polymorphism in this gene results in both coding and non-coding variants. The reference genome contains the non-coding allele. The non-functional allele of this gene is associated with elite athlete status.

**Antigen:** Recombinant protein of human ACTN3

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expressed only in a subset of type 2 skeletal muscle fibers.