

Anti-Human ARHGAP25 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA12922

3th Edition

Description: Rho GTPase activating protein 25 is a protein that in humans is encoded by the ARHGAP25 gene. The gene is also known as KAIA0053. ARHGAP25 belongs to a family of Rho GTPase-modulating proteins which are implicated in actin remodeling, cell polarity, and cell migration. GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.

Antigen: Recombinant protein of human ARHGAP25

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: