

Anti-Human/Mouse ARHGAP5 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA12924

3th Edition

Description: Rho GTPase activating protein 5 negatively regulates RHO GTPases, a family which Human, Mouse may mediate cytoskeleton changes by stimulating the hydrolysis of bound GTP. Two transcript variants encoding different isoforms have been found for this gene. GTPase-activating protein for Rho family members. May play a role in the reduction of the p21rasGTPase-activating potential of RASA1/p120GAP.

Antigen: Recombinant protein of human ARHGAP5

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in kidney, brain, liver and lung.