

Anti-Human/Mouse/Rat ANKRD54 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA12956

3th Edition

Description: Plays an important role in regulating intracellular signaling events associated with erythroid terminal differentiation. Interacts (via ankyrin repeat region) with LYN (via SH3-domain) in an activation-independent status of LYN By similarity. Forms a multiprotein complex with LYN and HCLS1 By similarity. Interacts with TSN2, VAV1, DBNL AND LASP1.

Antigen: Recombinant protein of human ANKRD54

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: