

**Anti-Human/Mouse ANGPTL7 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA12959**

---

3th Edition

**Description:**Angiopoietin-related protein 7 is a protein that in humans is encoded by the ANGPTL7 gene, a member of the ANGPTL family. Highly and specifically expressed in the cornea where is confined to the stromal layer.

**Antigen:**Recombinant protein of human ANGPTL7

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Highly and specifically expressed in the cornea where is confined to the stromal layer.