

**Anti-Human AGTRAP Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA12999**

---

3th Edition

**Description:** This gene encodes a transmembrane protein localized to the plasma membrane and perinuclear vesicular structures. The gene product interacts with the angiotensin II type I receptor and negatively regulates angiotensin II signaling. Alternative splicing of this gene generates multiple transcript variants encoding different isoforms.

**Antigen:** Recombinant protein of human AGTRAP

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Ubiquitous but more abundant in kidney, heart, pancreas and thyroid.