
CASP3, 176-277aa, Human, E.coli**Cat.NO.: TP01433**

3th Edition

Synonyms:Caspase 3, CPP32, CPP32B, SCA-1**Description:**CASP3 is a protein which is a member of the cysteine-aspartic acid protease (caspase) family. Sequential activation of caspases plays a central role in the execution-phase of cell apoptosis. Caspases exist as inactive proenzymes which undergo proteolytic processing at conserved aspartic residues to produce two subunits, large and small, that dimerize to form the active enzyme. This protein cleaves and activates caspases 6, 7 and 9, and the protein itself is processed by caspases 8, 9 and 10. It is the predominant caspase involved in the cleavage of amyloid-beta 4A precursor protein, which is associated with neuronal death in Alzheimer's disease. Recombinant human CASP3 protein, fused to His-tag at N-terminus, was expressed in E.coli.**Form:**Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.4M Urea, 10% glycerol**Molecular Weight:**12 kDa (103aa)**Sequences:**MSGVDDDMACHKIPVEADFLYAYSTAPGYYSWRNSKDGSWFIQSLCAMLKQYADKLEFMHILTRVNRKVATEFES
FSFDATFHAKKQIPCIVSMLTKELYFYH**Purity:**> 95% by HPLC**Concentration:**0.5 mg/ml (determined by Bradford assay)**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.