

**CD244, 19-224aa, Human, hlgG-His-tag, Baculovirus**

**Cat.NO.: TP01494**

3th Edition

**Synonyms:** Natural killer cell receptor 2B4 isoform 1, CD244, 2B4, NAIL, NKR2B4, Nmrk, SLAMF4

**Description:** CD244, also known as natural killer cell receptor 2B4 isoform 1, is a natural killer and T cell immunoglobulin superfamily surface protein and ligand of CD48. It is a cell surface glycoprotein related to CD2 and implicated in the regulation of natural killer and T lymphocyte function. It is expressed on all human NK cells, a subpopulation of T cells, basophils and monocytes. It activates NK cell mediated cytotoxicity, induces secretion of IFN-gamma and matrix metalloproteinases, and NK cell invasiveness. Recombinant human CD244, fused to hlgG-His-tag at C-terminus, was expressed in insect cell and purified by using conventional chromatography techniques.

**Form:** Liquid. In Phosphate Buffered Saline (pH 7.4) containing 10% glycerol.

**Molecular Weight:** 50.2kDa (448aa) 50-70kDa (SDS-PAGE under reducing conditions)

**Sequences:**

**Purity:** > 95% by HPLC

**Concentration:** 0.5mg/ml (determined by Absorbance at 280nm)

**Endotoxin Level:** <1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:** Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.