

CNTF, 1-200aa, Human, His-tag, E.coli (Bioactivity Validated)

Cat.NO.: TP01673

3th Edition

Synonyms:Ciliary neurotrophic factor

Description:CNTF (Ciliary neurotrophic factor) is a polypeptide hormone whose actions appear to be restricted to the nervous system where it promotes neurotransmitter synthesis and neurite outgrowth in certain neuronal populations. This protein is a potent survival factor for neurons and oligodendrocytes and may be relevant in reducing tissue destruction during inflammatory attacks. CNTF protein is structurally related to leukemia inhibitory factor (LIF), interleukin-6 (IL-6), interleukin-11 (IL-11) and oncostatin M (OSM). It is localized in the cell nucleus subsequent to receptor binding. Recombinant human CNTF, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

Form:Liquid. In 20 mM Tris-HCl buffer (pH 8.5) containing 1 mM DTT, 30% Glycerol, 0.2M NaCl.

Molecular Weight:25kDa (220aa), confirmed by MALDI-TOF

Sequences:

Purity:> 95% by HPLC

Concentration:1 mg/ml (determined by Bradford assay)

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.