

Recombinant Human IL-18 Receptor Accessory Protein/IL-18RAcP/IL-1R7 Protein(C-6His)

Cat.NO.: TP04617

3th Edition

Synonyms: Interleukin-18 receptor accessory protein (IL18RAP); Accessory protein-like; CD218 antigen-like family member B; CDw218b; IL-1R accessory protein-like; Interleukin-1 receptor 7; Interleukin-18 receptor accessory protein-like; Interleukin-18 receptor beta? CD218b and IL1R7.

Description: IL-18RAP is a single-pass type I membrane protein and contains two Ig-like C2-type domains and one TIR domain, IL18RAP can be induced by IFN-alpha and IL12 in nature killer cells and T-cells. The coexpression of IL18R1 and IL18RAP is required for the activation of NF-kappa B and JNK in response to IL-18. Defects in IL18RAP are associated with Coeliac disease.

Form: PBS

Molecular Weight: 39.4 kDa

Sequences: Phe20-Arg356

Purity: > 95% by HPLC

Concentration:

Endotoxin Level: <1.0 EU per 1 ug of protein (determined by LAL method)

Storage: Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.