

Recombinant Human TNF-Related Apoptosis-Inducing Ligand Protein

Cat.NO.: TP04708

3th Edition

Synonyms:Tumor Necrosis Factor Ligand Superfamily Member 10; Apo-2 Ligand; Apo-2L; TNF-Related Apoptosis-Inducing Ligand; Protein TRAIL; CD253; TNFSF10; APO2L; TRAIL

Description:Human TNFSF10 is a type II transmembrane protein with an intracellular N-terminus and a 'TNF homology domain' (THD) at the extracellular C terminus. TNFSF10 can interact with several distinct receptors. Two of these receptors that belongs to TNFR superfamily, DR4 (TRAIL-R1) and DR5 (TRAIL-R2/TRICK2), are plasma membrane proteins containing intracellular death domains essential for activating apoptosis. TNFSF10 is promising for cancer therapy because it is cytotoxic and activates apoptosis in the majority of malignant cells, but not in normal cells.

Form:PBS

Molecular Weight:19.5 kDa

Sequences:Arg115-Gly281

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.