

Recombinant Human Protein FAM3C Protein(C-6His)

Cat.NO.: TP04779

3th Edition

Synonyms:Protein FAM3C; Interleukin-Like EMT Inducer; FAM3C; ILEI

Description:FAM3C, also called interleukin-like EMT inducer, usually exist in most secretory epithelia. It belongs to the FAM3 family according to their sequence similarities. The up-regulation and/or mislocalization in breast cancer and liver carcinoma cells of FAM3C is strongly correlated with metastasis formation and survival. FAM3C can be involved in retinal laminar formation and promote epithelial to mesenchymal transition.

Form:PBS

Molecular Weight:23.2 kDa

Sequences:Gln25-Asp227

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.